



**sartorius**  
mechatronics

# Speed meets analytical precision



## Measure the moisture content in seconds with the new LMA200

Effective process monitoring demands high-speed analysis along with enhanced analytical precision.

Developed for fast at-line monitoring, the LMA200 uses high-energy microwave radiation to heat samples. The built-in analytical weighing system ensures that the measurements are performed with the required precision.

LMA200 – the fastest way to gravimetric moisture analysis!



[www.sartorius-mechatronics.com](http://www.sartorius-mechatronics.com)  
turning science into solutions